


ONE PERSON IN GROUP-> CLICK USE TEMPLATE



WHEN THE COPIED TEMPLATE OF THESE SLIDES OPENS IN A NEW TAB:

- ADD YOUR Group NAME TO THE TITLE
- SHARE WITH GROUP MEMBERS (USE SHARE BUTTON)
- COPY THE TITLE AND LINK TO YOUR RECORDS (Digital folder/table of contents)
- Remember it is advisable to make a PERSONAL copy and link it when DONE

Dilutions

ESSENTIAL QUESTIONS

Remember, you should be able to thoroughly answer these AFTER learning through following resources.

- How can we calculate dilutions and concentrations we need in the lab?
- How can we safely prepare dilutions for our needs?

Instructions

Follow through the **Engage to Apply** resources in order .

- Once complete, choose from Extend resources as interested.
- Record your learning in your **these slides- Remember to keep track of where to find these slides for studying!**
- Complete the **Essential Questions Response in the Apply Section** .
- Try the Extend section and optional Apply practice

Hyperdoc Elements

Click for quick access to:

- [Engage](#)
- [Explore](#)
- [Explain](#)
- [Apply](#)
- [Extend](#)

Terminology - As you encounter terms new to you, define in your own words here, include examples (make them up!)




Part 1: ENGAGE

Explore - Part A: 5 minutes

Make two observations & one inference about the image on the right:

Observations:

Inference:



Part 2: EXPLORE

Explore - Preparing Solutions

Preparing solutions with known concentrations is an important skill in chemistry

- Prepared Solutions with known concentrations - **Standard Solution**
- Purchased Solutions from a Chemical Provider - **Stock Solution**
 - Both used for accurate stoichiometric calculations

There are two ways to prepare Standard Solutions:


1. Dissolving solid solutes in a Volume of water
2. Diluting a Stock Solution to a required Concentration

Explore - Glassware & Tools


Because solutions need to be **accurate** and **reliable**, there are special pieces of glassware we can use to prepare solutions:

- 1) **Volumetric Flask:** a flask that measures a single volume with very high accuracy
- 2) **Graduated Cylinder:** a piece of glassware with gradations, that can measure many volumes
- 3) **Pipette:** a piece of glassware with gradations that can measure volumes with high accuracy


Explore - Glassware & Tools




Beaker: Used to contain and dispense liquids



Erlenmeyer flask: Used to contain and dispense liquids; shaped to slow evaporation, allows for a stopper, and lends itself to swirling liquids to mix them



Volumetric flask: Used to measure and contain a single, extremely accurate volume of liquid



Graduated cylinder: Used to measure and dispense liquid by pouring




Graduated pipette: Used to measure and dispense liquid using suction

Explore - Measuring Volumes

When measuring **volumes of liquids**, it is important to consider the **meniscus**

- ❖ When measuring liquids, **you must measure from the BOTTOM of the meniscus**



Explore - Preparing Solutions

Solutions with **Solid Solutes** can be prepared using the following steps:

- 1) Measure mass of solute using balance
- 2) Dissolve solute in water (less than total volume)
- 3) Pour mixture to **volumetric flask**
- 4) Rinse beaker with distilled water, and pour mixture into volumetric flask.
Repeat this step twice.
- 5) Add solvent to correct volume.
- 6) Mix and store solution in container.

Explore - Preparing Solutions



(a) An amount of solute is weighed out on an analytical balance and then

(b) A portion of the solvent is added to the volumetric flask.

(c) The mixture is swirled until *all* of the solute is dissolved.

(d) Additional solvent is added up to the mark on the volumetric flask.

Explore - Solutions from Dilutions

Highly concentrated solutions from a manufacturer are called stock solutions

- Concentrations often *too high* for experiments
- diluted to useful concentrations by adding more solvent

The **solution-dilution equation** can help us determine how to dilute to get desired concentrations

The diagram shows the solution-dilution equation $C_1 V_1 = C_2 V_2$ enclosed in a green rectangular box. Four red curved arrows point from text labels to the variables in the equation: 'Starting concentration' points to C_1 , 'Starting volume' points to V_1 , 'Final concentration' points to C_2 , and 'Final volume' points to V_2 .

Explore - Solution-Dilution Equation

The diagram illustrates the derivation of the dilution equation. On the left, two equations are shown, each preceded by a large blue chevron pointing right:

$$C_1 = \frac{n_1}{V_1} \quad \text{and} \quad C_2 = \frac{n_2}{V_2}$$

These lead to the equations:

$$n_1 = V_1 \times C_1 \quad \text{and} \quad n_2 = V_2 \times C_2$$

A large blue curly brace in the center groups these two equations, with the text $n_1 = n_2$ written inside it. To the right of the brace is the dilution equation:

$$V_1 C_1 = V_2 C_2$$

Four blue arrows point from text labels to the terms in the dilution equation:

- An arrow from "Concentration of standard" points to C_2 .
- An arrow from "Volume of standard" points to V_2 .
- An arrow from "Concentration of stock" points to C_1 .
- An arrow from "Volume of stock" points to V_1 .

Discuss! Why does $n_1 = n_2$? Remember you are DILUTING a volume (v_1) of the stock to a new higher volume but lower concentration solution.

Explore - Solution-Dilution Practice

Ms.K & Dr. Loney want to prepare 2.00L of a 0.100M solution of HCl. The concentration of the stock 12.0M. What volume of stock HCl is needed to make this solution?


$$C_1V_1 = C_2V_2$$

$$V_1 = \frac{C_2V_2}{C_1} = \frac{(0.100M)(2.00L)}{12.0M} = 0.0167 \text{ L} = 16.7\text{mL}$$

Explore - Labelling Solution Containers

Solutions should have the following information on their containers when stored:


- Chemical name of solute
- Concentration
- Date prepared
- Name of person who prepared the solution



Explore - Dissolving



{Optional} For further help you can watch this video



Part 3: Apply

Apply- Complete the questions



Show all your work! Include formulas, conversions, units, etc.

- 1) For a lab, it requires 305mL of a 4M sugar solution. The lab then tells you to add 695mL of water to the solution. What is the final concentration of the sugar solution?
- 2) You want to create 375mL of a 0.450M solution of Sulfuric Acid. If you have a stock solution of 12M, how many mL of the stock is required?
- 3) 280mL of water is added to 50mL of a 2M solution of Lead Chloride. What is the new concentration of the solution?

Next few slides provide answers to activities

Apply- Complete the questions



Show all your work! Include formulas, conversions, units, etc.

- 1) For a lab, it requires 305mL of a 4M sugar solution. The lab then tells you to add 695mL of water to the solution. What is the final concentration of the sugar solution? **Ans: 1.22M**
- 2) You want to create 375mL of a 0.450M solution of Sulfuric Acid. If you have a stock solution of 12M, how many mL of the stock is required? **Ans: 14mL**
- 3) 280mL of water is added to 50mL of a 2M solution of Lead Chloride. What is the new concentration of the solution? **Ans: 0.3M**

Part 5: Beyond the course extensions {optional}